



- ▶ Mobile disinfection station
- ✓ Disinfectant dispenser device



- ▶ Available as an option
- ✓ Bartscher disinfectant dispenser PS 0.9L-W
- ✓ Bartscher drip tray PS 0.9L-W



The disinfection station is the ideal mobile solution for hand hygiene, particularly in entrance areas. It does not need to be mounted to a wall or to the floor. With the optional disinfectant dispenser and drip tray, the disinfection station enables top hygiene standards for pubs, restaurants, hotels, and petrol stations etc.

- Designed for: Bartscher Disinfectant dispenser PS 0.9L-W
Bartscher Drip tray PS 0,9L-W
- Disinfectant dispenser device: Yes
- Delivery state: Building kit
- Properties: -
- Not included in delivery: Disinfectant dispenser
Drip tray
- Important information: -
- Including: Attachment material
- Material: Stainless steel
- Glove dispenser device: No

- Integrated drip tray: No
- Model: Upright unit
- Colour: Silver
- Set comprises: 1 plinth
1 column
- Size: W 300 x D 400 x H 1,170 mm
- Weight: 7.4 kg

Disinfection station D1 1170

Add on Products

Disinfectant dispenser PS 0.9L-W

0,9 L




- Content: 0,9 litres
- Can be refilled: Yes
- Filling level check: Transparent plastic container
- Operation: One-hand operation
Elbow operation
Manual
- Lockable: No
- Properties: Cartridge-free
- Including: 1 wall holder
Attachment material
- Not included in delivery: Filling

Code-No. 850019
GTIN 4015613704746

Drip tray PS 0,9L-W

STAINLESS
STEEL



- Material: Stainless steel
Polypropylene
- Important information: -
- Properties: Collection tray suitable for dishwasher
- Collection tray: Size: W 147 x D 240 x H 40 mm
- Size: W 147 x D 285 x H 377 mm
- Weight: 0.46 kg

Code-No. 850018
GTIN 4015613699431

Disinfectant LEOCID® P7



- Order quantity unit: 1 box (12 bottles)
- Content: 12 x 1 litre
- Model: Liquid
- Ambient temperature: 5 °C to 25 °C
- Ph level: 7
- Properties: Partially viricidal
Bactericidal, yeasticidal, tuberculocidal
Perfume-free
- Designed for: Disinfectant dispenser
- Application: Allow 3 ml of undiluted liquid to act for 30 seconds (as recommended by RKI)

Code-No. 173084
GTIN 4029888031308